## **Glucose CSF**

Synonyms	CSF Glucose, Spinal fluid glucose
Abbreviation	CSF GLU
Specimen type	Sterile CSF
Container	CSF tube
Volume	1 mL Minimum: 0.5 mL
Collection requirements	Deliver to Lab immediately; must arrive within 6 hours. Refrigerate if analysis will be delayed.
S/H instructions	Refrigerate.
Stability	3 days at room temp. 1 month refrigerated.
Reference range	40 – 70 mg/dL
Routine / STAT T.A.T	Within 4 hours / 1 hour
Section	Chemistry
LAB code	01470