

Troponin High Sensitive

Synonyms	
Abbreviation	TROP hsTrop
Specimen type	Plasma Serum can be used, but specimen type must be the same for all tests within a series.
Container	Green/Yellow (PST) , Green/White or Dark Green Vacutainer.
Volume	1 mL plasma Minimum: 0.5 mL Plasma
Collection requirements	Green Top is preferred sample.
S/H instructions	
Stability	Test within 8 hours. Centrifuged: 72 hours refrigerated (aliquot). 30 days frozen (aliquot); thaw only once.
Reference range	Male \leq 35 pg/mL Female <15 pg/mL
Routine / STAT T.A.T	Within 4 hours / 1 hour
Section	CHEM
LAB code	08275