

Type + Screen Panel

Synonyms	ABO/RH + Screen, T and S
Abbreviation	TANDS
Specimen type	Whole blood Preferred: Plasma & cells from Pink top Acceptable: Plasma & cells from Lavender
Container	2 6mL Pink EDTA (Collect 4 pink tubes if patient has an RBC antibody history) Acceptable: Lavender
Volume	6 mL Minimum: 3 mL
Collection requirements	Avoid hemolysis
S/H instructions	Refrigerate
Stability	7 days at 1-8 C
Reference range	NA
Routine / STAT T.A.T	Within 4 hours / 1 hour
Section	Blood Bank
LAB code	03029