

URIC ACID

Synonyms	URAC
Abbreviation	URAC
Specimen type	Preferred: Plasma Acceptable: Serum
Container	Green/Yellow (PST), Dark Green, Green/White Vacutainer, or Green Microtainer; Red/Yellow (SST) or Plain Red Vacutainer or Red Microtainer.
Volume	1 mL Minimum: 0.5 mL Plasma/Serum
Collection requirements	
S/H instructions	
Stability	Centrifuged: 3 days at room temperature. 7 days refrigerated. 6 months frozen.

Reference range	<table border="1"> <thead> <tr> <th>AGE</th> <th>RANGE</th> </tr> </thead> <tbody> <tr> <td>Child</td> <td>2.0-5.5 mg/dL</td> </tr> <tr> <td>Adult Male</td> <td>3.9-7.8 mg/dL</td> </tr> <tr> <td>Adult Female</td> <td>2.2-6.8 mg/dL</td> </tr> </tbody> </table>		AGE	RANGE	Child	2.0-5.5 mg/dL	Adult Male	3.9-7.8 mg/dL	Adult Female	2.2-6.8 mg/dL
	AGE	RANGE								
	Child	2.0-5.5 mg/dL								
	Adult Male	3.9-7.8 mg/dL								
Adult Female	2.2-6.8 mg/dL									
Routine / STAT T.A.T	Within 4 hours / 1 hour									
Section	CHEM									
LAB code	01700									