

## Sars-CoV-2 (Covid 19) + Influenza A&B + RSV by PCR

<b>Synonyms</b>	RSV, Flu, Influenza, Covid 19
<b>Abbreviation</b>	NA
<b>Specimen type</b>	Preferred: Nasopharyngeal Swab collected In universal transport media. Alternate Sources: Nasal
<b>Container</b>	Universal or Viral transport media
<b>Volume</b>	NA
<b>Collection requirements</b>	Not to be confused with Aspirin Platelet Function Assay.
<b>S/H instructions</b>	
<b>Stability</b>	Room Temp: 48 hrs Refrigerated (2-8C): 7 Days
<b>Reference range</b>	Negative
<b>Routine / STAT T.A.T</b>	Within 4 hours / 1 hour
<b>Lab Code</b>	13969